

Malleo Sprint

The number one choice for chronic instability



Quality for life



Malleo Sprint

The new Malleo Sprint ankle orthosis is the number one therapy choice in case of chronic ligament instability in the foot/ankle region. Thanks to the special plastic exoskeleton and strapping, the Malleo Sprint assures a high degree of stabilisation without having a negative impact on plantar and dorsal flexion. Undesirable eversion and inversion of the foot are limited.

When it is closed, the orthosis structure, which adapts perfectly to the contours of the body, provides targeted, stable support and limitation of movement. Breathable materials and excellent padding ensure a high level of wearer comfort.

Features and Benefits

A high level of wearer comfort thanks to an outstanding anatomical fit

Anatomically-shaped padding assures a high level of wearer comfort, while graduated padding underneath the foot prevents pressure points.

A high degree of stability

The combination of a plastic skeleton and strapping provides a high level of stability and prevents unwanted movements.

Excellent Support

The integrated hook and loop system optimises the support provided by the laces and allows for easy donning and doffing.

Straightforward application

An additional loop further supports donning and doffing.


Indications

- Chronic ankle joint instability
- Functional treatment of ligament ruptures
- Severe ankle distortions



► The orthosis is universal for the left or right side.

Order no.	Size	Circumference	
		(cm)	(inch)
50S3=XS-1	XS	25,5–28,0	10–11
50S3=S-1	S	28,0–30,5	11–12
50S3=M-1	M	30,5–33,0	12–13
50S3=L-1	L	33,0–35,5	13–14
50S3=XL-1	XL	35,5–38,0	14–15



Otto Bock HealthCare GmbH
Max-Näder-Straße 15 · 37115 Duderstadt
T +49 5527 848-1424 · F +49 5527 848-72330
export@ottobock.de · www.ottobock.com